

## Silica Gel blue 2-4 mm with indicator, beads

[Product in eShop](#)

[Home](#) / [Brands](#) / [CHEMSOLUTE](#) / Silica Gel blue 2-4 mm with indicator, beads

### Product Data

Drying agent

Alternative product: environmentally friendly silica gel.

- Patented, organic humidity indicator - **heavy metal free**
- Reversible colour change

#### Note:

- We recommend activating silica gels before using it for the first time

#### Properties:

- Free of cobalt(II) chloride
- Regeneration in drying cabinet at approx. 80-100 °C for approx. 2-4 hours
- Colour change from blue to fuchsia

<b>Formula</b>	O <sub>2</sub> Si
<b>CAS No.</b>	7631-86-9
<b>EC No</b>	231-545-4
<b>LGK</b>	13

Item No.	Content	Packaging	Price
2440.2500	2.5 Kilo		CHF 125.30

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

[www.huberlab.ch](http://www.huberlab.ch)