

# Phenol equilibrated, stabilized

 Product in eShop

Home / [Assortment](#) / [Chemicals](#) / [Biochemicals](#) / [Further](#) / Phenol equilibrated, stabilized

## Product Data

<b>CAS No.</b>	108-95-2
<b>HS-No</b>	29071100
<b>EC No</b>	203-632-7
<b>Lagerung</b>	2-8°C
<b>LGK</b>	6.1A
<b>Hazard Statement(s)</b>	H301+H311+H331-H314-H341-H373
<b>Safety Instructions</b>	P280-P301+P330+P331- P302+P352-P305+P351+P338- P309-P310
<b>Signal Word</b>	Gefahr
<b>Class / PG</b>	6.1/II
<b>UN-No.</b>	UN2821
<b>WHC</b>	2



Item No.	Content	Packaging	Price
A1153.0500	500 milliliter	Glasflasche	CHF 187.45

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 17.04.2026 at 21:27 o'clock.