

Thermo shaker uniTHERMIX pro LLG- Labware

[Product in eShop](#)

Home / [Labortops_Q2_2025](#) / Thermo shaker uniTHERMIX pro LLG-Labware

Product Data

The LLG-uniTHERMIX pro is ideal for shaking and temperature control of a variety of tubes or plates (depending on the block used). With its high temperature uniformity, fast heating rate and excellent mixing results the LLG-uniTHERMIX pro is perfect for the incubation of samples of any kind, in particular for Life Science applications, e.g. Enzyme reactions, transformations, denaturation of DNA and RNA. The heating blocks are magnetically attached to the LLG-uniTHERMIX pro and do not have to be screwed down.

- Easy and quick change of heating blocks due to magnetic fastening
- Precise temperature control system
- Temperatures up to 100 °C
- Temperature accuracy of ± 0.5 °C (20 ... 45 °C)
- Efficient mixing up to 1500 rpm (depending on the heating block)
- Up to 9 storable programs
- Many block options available
- Automatic block recognition

LLG-uniTHERMIX 1 pro: with shaking and heating function, Temperature range: RT + 5 °C ... approx. 100 °C

LLG-uniTHERMIX 2 pro: with shaking, cooling and heating function, Temperature range: 15 °C below RT ... approx. 100 °C

Blocks must be ordered separately.

Scope of supply: thermal shaker, 4 suction feet

Item No.	Description	UNSPSC	ECLASS	Price
600.2634.82	LLG-uniTHERMIX 1 pro Thermomixer with heating, with EU-plug		32030000	CHF 1'500.00

600.2634.84	LLG-uniTHERMIX 2 pro Thermomixer with heating and cooling, with EU- plug	32030000	CHF 1'965.00
600.2634.87	Heating block for LLG-uniTHERMIX 1/2 pro used for 1,5ml tubes, 24 holes	32030000	CHF 225.00
600.2634.94	Heating block for LLG-uniTHERMIX 1/2 pro used for 2ml tubes, 24 holes	32030000	CHF 225.00
600.2634.95	Heating block for LLG-uniTHERMIX 1/2 pro used for 5ml tubes, 8 holes	32030000	CHF 222.70
600.2634.96	Heating block for LLG-uniTHERMIX 1/2 pro used for 15ml tubes, 8 holes	32030000	CHF 380.00
600.2634.97	Heating block for LLG-uniTHERMIX 1/2 pro used for 50ml tubes, 4 holes	32030000	CHF 380.00
600.2634.98	Heating block for LLG-uniTHERMIX 1/2 pro, used for 0.2 ml/96 well PCR plate	32030000	CHF 221.80
600.2634.99	Heating block for LLG-uniTHERMIX 1/2 pro, used for 96/384 microplate	32030000	CHF 225.00

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 02.05.2026 at 08:00 o'clock.