

Glove dispenser, acrylic glass LIG-Labware



Product in eShop



Home / [Brands](#) / [LLG-Labware](#) / Glove dispenser, acrylic glass LIG-Labware

Product Data

Glove dispensers are suitable for glove packs with the dimensions approx. 240 x 75 x 125 mm, with front insert or side slot. High-quality, sturdy design. Mounting material included.

Item No.	Description	UNSPSC	ECLASS	Price
600.4053.72	LLG-Glove dispenser for 1 box side filling acrylic glass		32030000	CHF 17.70
600.4053.73	LLG-Glove dispenser for 3 boxes front insert, acrylic glass, with wall mounting material		32030000	CHF 50.00
600.4053.74	LLG-Glove dispenser for 3 boxes lateral insert, acrylic glass		32030000	CHF 52.50

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 02:59 oclock.

www.huberlab.ch