

# Magnetic stirring bar retriever PTFE, Ø 8 mm LLG-Labware



Product in eShop

Home / [Labortops\\_Q3\\_2025](#) / Magnetic stirring bar retriever PTFE, Ø 8 mm LLG-Labware

## Product Data

PTFE-encapsulated stirring bar retriever with strong permanent magnet (Alnico 5), universal chemical resistance. For the removal of stirring bars from aggressive liquids, prevents loss of stirring bars. Lower end diameter 10 mm.

Item No.	Description	UNSPSC	ECLASS	Price
600.0085.01	LLG-Magnetic stirring bar retriever 270 mm PTFE, Ø 8 mm		32030000	CHF 13.00 <del>CHF 21.10</del>
600.1977.17	LLG-Magnetic stirring bar retriever 150 mm PTFE, Ø 8 mm		32030000	CHF 10.30 <del>CHF 16.70</del>
600.1977.25	LLG-Magnetic stirring bar retriever 250 mm PTFE, Ø 8 mm		32030000	CHF 10.95 <del>CHF 17.80</del>
600.1977.35	LLG-Magnetic stirring bar retriever 350 mm PTFE, Ø 8 mm		32030000	CHF 11.85 <del>CHF 19.25</del>

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 06.05.2026 at 22:34 o'clock.