

## Pipette controller easy LLG-Labware



Product in eShop

Home / [Labortops Q2 2025](#) / Pipette controller easy LLG-Labware

### Product Data

The LLG-Pipette controller *easy* enables fast and exact pipetting with glass and plastic pipettes.

- Ergonomic and lightweight construction for constant ease of use over prolonged usage
- Unique nose cone angle for enhanced flexibility in operation e.g. in fume hoods
- Eco-friendly rechargeable NiMH batteries, offering up to 8 hours of continuous working
- Easily replaceable batteries
- Autoclavable tip cone for avoiding cross-contamination

**Scope of supply:** Pipette controller, manual, charger, wall-mount, desktop stand

Item No.	type no.	Price
600.2636.83	Pipette controller easy	CHF 169.50

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 06.05.2026 at 21:51 o'clock.