


## Slide mailer, high density PP for slides LLG-Labware

[Product in eShop](#)

Home / [Labortops\\_Q1\\_2025](#) / Slide mailer, high density PP for slides LLG-Labware

### Product Data

Manufactured from high density polypropylene, for shipping, storing or staining slides. Flat mailers can hold 1, 2, 3 or 5 standard (up to 76mm x 26mm) slides securely.

- Light weight and flat
- With marking area and "Top" printed on top surface
- With large, easy-to secure snaps.
- Slotted channels inside mailer keep slides safely separated

Item No.	Description	UNSPSC	ECLASS	Price
600.1612.91	LLG-Slide mailer, high density PP for 1 slide, pack of 10		32030000	CHF 11.50
600.1612.92	LLG-Slide mailer, high density PP for 2 slides, pack of 10		32030000	CHF 13.15
600.1612.93	LLG-Slide mailer, high density PP for 3 slides, pack of 10		32030000	CHF 14.50
600.1612.94	LLG-Slide mailer, high density PP for 5 slides, pack of 10		32030000	CHF 14.70

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 03.05.2026 at 03:00 o'clock.