

Pipette tips Macrotips Qualitix® Socorex

[Product in eShop](#)[Home](#) / [Labortops_Q3_2025](#) / [Pipette tips Macrotips Qualitix® Socorex](#)

Product Data

- Superior pipette tips cover a range between 0.2 mL and 2 mL/5 mL/10 mL.
- All models autoclavable.
- Optimal tightness fitting on Socorex® pipettes and compatibility with other brands.

Item No.	description	Price
2.7310.12	Autoclaving box with sliding lid, empty, for 2 mL macrotips	CHF 45.50
2.7312.02	macrotips	CHF 94.20
2.7312.03	---	CHF 1'131.40
2.7312.06	macrotips, graded	CHF 100.00
2.7312.07	Macrotip, graded	CHF 71.70
2.7312.08	---	CHF 1'061.50
2.7312.10	macrotips	CHF 106.50
2.7312.11	Empty rack for 24x 10 mL macrotips	CHF 69.00
2.7312.12	---	CHF 1'134.60

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 19.04.2026 at 07:09 o'clock.