

Samplepack Headspace Agilent

[Product in eShop](#)[Home](#) / [Brands](#) / [CG](#) / [Samplepack Headspace Agilent](#)

Product Data

- 8.1802.34: 20 ml Precision Thread Vial, 75.5 x 22.5 mm, clear glass, 1st hydrol.class, rounded bottom
- 8.1801.77: 18 mm UltraClean Seal: magnetic screw cap, with hole, Silicone white/PTFE blue, 55° shore A, 1.5 mm
- 8.2016.72: 20 ml Headspace-Vial, 75.5 x 22.75 mm, clear glass, 1st hydrol. class, with bevelled HS-crimp neck, flat bottom
- 8.2002.36: 20 ml Headspace-Vial, 75.5 x 22.5 mm, clear glass, 1st hydrol.class, DIN-crimp neck, long neck, flat bottom
- 8.2006.40: 20 mm UltraClean aluminium crimp cap, clear lacquered, with hole, Silicone white/PTFE beige, 45° shore A, 3 mm
- 8.2001.33: 20 mm UltraClean aluminium crimp cap, clear lacquered, with hole, Silicon blue transparent/PTFE white, 45° shore A, 3 mm

Item No.	Price
8.2018.07S	for free

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 01.05.2026 at 11:25 o'clock.