

Samplepack Headspace Shimadzu

[Product in eShop](#)[Home](#) / [Brands](#) / [CG](#) / [Samplepack Headspace Shimadzu](#)

Product Data

- 8.1802.34: 20 ml Headspace Precision Thread Vial, 75.5 x 22.5 mm, clear glass, 1st hydrol. class, rounded bottom
- 8.1805.75: 18 mm UltraClean Magnetic Universal Screw Cap, silver, with hole, Silicone blue transparent/PTFE white, 45° shore A , 1.3 mm
- 8.1805.77: 18 mm UltraClean Magnetic Universal Screw Cap, silver, with hole, Silicone white/PTFE blue, 55° shore A , 1.5 mm
- 8.2002.37: 20 ml Headspace-Vial, 75.5 x 22.5 mm, Klarglas, 1st hydrol. class, with DIN-crimp neck, rounded bottom
- 8.2008.36: 20 mm UltraClean Magnetic Crimp Cap, silver, with hole, Silicone blue transparent/ PTFE white, 45° shore, 3 mm
- 8.2002.15: 20 mm UltraClean Magnetic Bimetall-Crimp Cap, blue, with hole, Silicone blue transparent/PTFE white, 45° shore, 3 mm

Item No.	Price
8.2018.10S	for free

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 03.05.2026 at 14:09 o'clock.