

## b.safe Tubing Holder Bohlender



Product in eShop

Home / [Brands](#) / [b.safe](#) / b.safe Tubing Holder Bohlender

### Product Data

For fixation of capillaries and tubing to avoid disarranging or bending of tubes. Self-adhesive bottom side. A b.safe Tubing Holder can arrange at once up to 8x O. D. 1,6 mm / 3x O. D. 3,2 mm or 1x O. D. 5 mm. Packing unit: pack with 10 pieces.

### Applications

To lead and fix tubing and capillaries without being disarranged or bent.

Item No.	length	height	width	Price
94.M562.01	30 mm	10 mm	11 mm	CHF 20.00

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 05.05.2026 at 07:46 o'clock.