

## b.safe Laboratory Flasks Bohlender

[Product in eShop](#)

Home / [Brands](#) / [b.safe](#) / b.safe Laboratory Flasks Bohlender

### Product Data

Laboratory flask made of Borosilicate glass including screw cap and spout made of PP.  
Packing unit: 1 piece.

Item No.	volume	thread	Price
94.R100.45	1000 ml	GL45	CHF 28.60
94.R105.80	1000 ml	GLS80	CHF 37.00

### More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 21.04.2026 at 04:55 o'clock.