

## Portable glove bag - mobile isolation enclosure Erlab

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Safety & Hygiene](#) / Portable glove bag - mobile isolation enclosure Erlab

### Product Data

#### **Erlab portable glove bag. A mobile isolation enclosure for working anywhere.**

Whether in the lab or in the field, set up in just a few minutes to carry out your research. Made from high-quality transparent PVC, the Erlab isolation enclosure is a multipurpose, mobile disposable “glove bag,” designed to protect people as well as the research and transport of biological materials. Highly portable, lightweight, and assembled in seconds, it allows you to work anywhere in complete protection. Its inclined shape creates an ergonomic working position. Its optimized seal makes it ideal for any field-based research, as well as for highly secure lab work. Inflated using compressed air at a pressure of 2.5 mm Hg, all enclosures are extremely resistant and are tested before being offered for sale. They come with a certificate of compliance and each enclosure can be identified by a unique serial number.

- In form of a pyramid
- Dimensions: W: 860 mm / D: 560 mm / H: 725 mm

Item No.	Price
17.0391.22	CHF 534.80

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 23.04.2026 at 08:28 o'clock.