

Disposable Reagent reservoirs economy, PS LLG-Labware

[Product in eShop](#)

Home / [Labortops_Q3_2025](#) / Disposable Reagent reservoirs economy, PS LLG-Labware

Product Data

Made from polystyrol. Designed for use with multi-channel pipettes and can accommodate up to 12-channel standard pipettes. Angled sides and a V-shaped bottom allow maximum recovery of reagents. All reservoirs have internal graduated volume lines and pour-spouts on the corners. 25 ml or 50 ml versions. Available as sterile and non-sterile versions.

- DNase-, RNase- and pyrogen-free certified

Item No.	packaging	volume	Price
600.2821.44	non-sterile, bulk	25 ml	CHF 28.10
600.2821.45	sterile, single packed	25 ml	CHF 53.10
600.2821.46	non-sterile, bulk	50 ml	CHF 36.85
600.2821.47	sterile, single packed	50 ml	CHF 75.05

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 22.04.2026 at 01:40 o'clock.