

## Agarose Standard LLG-Labware

[Product in eShop](#)

[Home](#) / [Labortops\\_Q1\\_2025](#) / Agarose Standard LLG-Labware

### Product Data

LLG-Agarose is ideal for routine separation of DNA and RNA fragments, as well as PCR products with electrophoresis. High transparency and low background of the agarose gel ensures clear and sharply defined bands.

- Very good dissolution and rapid gel formation
- High gel strength even at low agarose concentrations
- Suitable for agarose concentrations from 0.75 up to 2 %
- Good separation efficiency for nucleic acid fragments between 250 bp and 23 kb
- Free of DNase and RNase

Appearance: White powder

Clarity (NTU): 4 (1.5 % gel)

Gelling point:  $36 \pm 1.5$  °C (1.5 % gel)

Melting point:  $88 \pm 1.5$  °C (1.5 % gel)

Gel strength: 1000 g/cm<sup>2</sup> (1 % gel)

Sulphate: 0.15 %

Item No.	Description	UNSPSC	ECLASS	Price
600.6815.11	LLG-Agarose Standard, 100g for DNA/RNA electrophoresis	41100000	32030000	CHF 66.60
600.6815.12	LLG-Agarose Standard, 500g for DNA/RNA electrophoresis	41100000	32030000	CHF 305.00

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 19.04.2026 at 12:32 oclock.

[www.huberlab.ch](http://www.huberlab.ch)