

Micropestles, stainless steel LLG-Labware

[Product in eShop](#)

[Home](#) / [Brands](#) / [LLG-Labware](#) / [Micropestles, stainless steel LLG-Labware](#)

Product Data

Polished stainless steel, micropestles. Precisely fitted for Eppendorf Safe-Lock Tubes™, manufactured in one-piece, smooth handle sections, can be disinfected easily. Minimum homogenate displacement due to tapered piston.

For Reaction tubes:

For homogenisation of cells, tissues or foodstuffs / feed samples.

- For tubes with 0.5 ml, 1.5 ml, and 2.0 ml
- Can be cooled down to -196 °C
- Sterilisable in an autoclave or dry heat at 200 °C
- Suitable for PCR applications
- Alternative to single-use plastic pestles
- Perfectly suitable for Eppendorf Safe-Lock Tubes

Item No.	Description	UNSPSC	ECLASS	Price
600.3145.00	LLG-micro-pestle stainless steel, for 0,5 ml vials	41100000	32030000	CHF 29.55
600.3145.01	LLG-micro-pestle stainless steel, for 1,5 ml vials	41100000	32030000	CHF 31.85
600.3145.02	LLG-micro-pestle stainless steel, for 2,0 ml vials	41100000	32030000	CHF 34.10

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 02.05.2026 at 08:17 o'clock.

www.huberlab.ch