

## Timer pro LLG-Labware



Product in eShop

Home / [Labortops\\_Q3\\_2025](#) / Timer pro LLG-Labware

### Product Data

**The LLG-Timer pro can be used simultaneously as real time clock and timer.** The countdown timer also runs when time is visible.

- Big 6-digit LCD display
- Digital stopwatch (counting up)
- Short timer (counting down)
- With memory function
- Stopwatch and countdown up to 23 h 59 min. 59 sec.
- Optical and acoustic alarm when timer expires
- Incl. magnet for mounting, hanging bracket and stand (fold-out)
  
- Dimensions (W x H x D): 82 x 68 x 20,5 mm
- Weight (with battery): 82 g
- Power supply: 2 x AAA 1,5 V (not included)
- Warranty: 12 months
- Measuring capacity (hr:min:sec) : 23:59:59

Item No.	type no.	Price
600.3113.51	LLG-Timer pro	CHF 21.30

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 28.04.2026 at 04:18 o'clock.