

Measuring cylinders, borosilicate glass 3.3, tall form, class A LLG- Labware

[Product in eShop](#)

Home / [Labortops Q1 2025](#) / Measuring cylinders, borosilicate glass 3.3, tall form, class A LLG-Labware

Product Data

Tall form. Conformity certified with batch ID.

- DIN EN ISO 4788
- With hexagonal base and spout
- Blue graduation

Item No.	for nominal volume	division	tol. +/-	Price
600.6862.02	5 ml	0.1 ml	±0.05 ml	CHF 19.60
600.6862.03	10 ml	0.2 ml	±0.10 ml	CHF 19.60
600.6862.04	25 ml	0.5 ml	±0.25 ml	CHF 20.80
600.6862.05	50 ml	1 ml	±0.50 ml	CHF 22.20
600.6862.06	100 ml	1 ml	±0.50 ml	CHF 26.50
600.6862.07	250 ml	2 ml	±1.00 ml	CHF 42.50
600.6862.08	500 ml	5 ml	±2.50 ml	CHF 58.80
600.6862.09	1000 ml	10 ml	±5.00 ml	CHF 38.10
600.6862.10	2000 ml	20 ml	±10.00 ml	CHF 64.70

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 11:25 o'clock.