


Volumetric flasks, borosilicate glass 3.3, class A LLG-Labware

[Product in eShop](#)

Home / [Labortops_Q2_2025](#) / Volumetric flasks, borosilicate glass 3.3, class A LLG-Labware

Product Data

Conformity certified with batch ID.

- DIN EN ISO 1042
- With standard ground joint and PE stopper
- Blue ring mark and inscription

Item No.	for nominal volume	Ground Joint	tol. +/-	Price
600.6862.32	5 ml	NS10/19	±0.040 ml	CHF 21.80
600.6862.33	10 ml	NS10/19	±0.040 ml	CHF 21.20
600.6862.34	20 ml	NS10/19	±0.040 ml	CHF 22.00
600.6862.35	25 ml	NS10/19	±0.040 ml	CHF 22.10
600.6862.36	25 ml	NS12/21	±0.060 ml	CHF 24.20
600.6862.37	50 ml	NS14/23	±0.100 ml	CHF 27.60
600.6862.38	50 ml	NS12/21	±0.060 ml	CHF 24.70
600.6862.39	100 ml	NS12/21	±0.100 ml	CHF 30.00
600.6862.40	100 ml	NS14/23	±0.100 ml	CHF 28.00
600.6862.41	200 ml	NS14/23	±0.150 ml	CHF 35.40
600.6862.42	250 ml	NS14/23	±0.150 ml	CHF 35.90
600.6862.43	500 ml	NS19/26	±0.250 ml	CHF 47.60
600.6862.44	1000 ml	NS24/29	±0.400 ml	CHF 36.00
600.6862.45	2000 ml	NS29/32	±0.600 ml	CHF 50.40

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

www.huberlab.ch