

## ESD clean room chair Werk Sitz



Product in eShop

Home / [Labortops Q1 2025](#) / ESD clean room chair Werk Sitz

### Product Data

- Backrest upholstered with leatherette, colour: black, (W top/bottom x H: 340/445 x 515 mm), backrest tilt: -10° to +8°
- Seat upholstered with leatherette, colour: black, (W back/front x D: 480/510 x 470 mm), seat tilt: -4° to +10°
- Column, chrome-plated, with column cover
- Five star base, polished with conductive dual castors
- Conductive chair for discharging static electricity electrical resistance:  $6 \times 10^6$  Ohm
- Certified for class-5 clean rooms acc. to DIN EN ISO 14644-1
- Load capacity 120 kg
- 5 years warranty

Item No.	Description	UNSPSC	ECLASS	Price
600.6615.13	Cleanroom chair WS 1720 RR ESD imitation leather, black, height adjustable cleanroom	41100000	32030000	CHF 1'035.00

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 24.04.2026 at 07:35 oclock.

[www.huberlab.ch](http://www.huberlab.ch)