

Autoclavable Bags, PP, with Biohazard printing LLG-Labware

[Product in eShop](#)

[Home](#) / [Labortops_Q1_2025](#) / Autoclavable Bags, PP, with Biohazard printing LLG-Labware

Product Data

Autoclavable bags imprinted with biohazard labelling. Made of polypropylene. Strong side seams for securely holding items. Front write-on panel. Autoclavable.

Item No.	w x l	colour	autoclavable	Price
600.4040.50	310 x 660 mm	transparent	at 121 °C	CHF 30.00
600.4040.51	415 x 600 mm	transparent	at 121 °C	CHF 33.00
600.4040.52	610 x 810 mm	transparent	at 121 °C	CHF 66.50

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

www.huberlab.ch