

Microcapillaries Calibre Scientific

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Cell Manipulation](#) / [Microcapillaries Calibre Scientific](#)

Product Data

It offers no magic bullet, but the right microcapillary for your cell manipulation application. Our high quality capillaries harmoniously complement our cell manipulation instruments, thereby providing you with the best set of tools for successful experimentation.

Product Information

A wide range of excellent microcapillaries, designed to give you fast, efficient and reproducible results for the most common applications. All microcapillaries offer you reproducible quality through narrowly defined specifications and intensive quality control, as well as the greatest security through effective sterilization methods.

Features

- TransferTips for stem cell transfer etc.
- Femtotips and Femtotip II for microinjection of minimal volumina
- Piezo-Drill Tips for piezo-assisted micromanipulation





Item No.	type no.	description	Price
5195000001	TransferTip® F (ICSI)	for sperm injection using the ICSI technique (for research use only), 35 ° tip angle, 4 µm inner diameter, 0.4 mm flange, sterile, 25 pcs.	CHF 721.00
5195000010	TransferTip® RP (ICSI)	for sperm injection using the ICSI technique (for research use only), 35 ° tip angle, 4 µm inner diameter, 0.5 mm flange, sterile, 25 pcs.	CHF 720.00
5195000028	TransferTip® R (ICSI)	for sperm injection using the ICSI technique (for research use only), 35 ° tip angle, 4 µm inner diameter, 1 mm flange, sterile, 25 pcs.	CHF 720.00

5195000036	VacuTip I	holding capillary (for research use only), 35 ° tip angle, 15 µm inner diameter, 100 µm outer diameter, 1 mm flange, sterile, 25 pcs.	CHF 549.00
5195000044	VacuTip II	holding capillary (for research use only), 25 ° tip angle, 60 µm inner diameter, 110 µm outer diameter, 0.5 mm flange, sterile, 25 pcs.	CHF 562.00
5195000079	TransferTip® (ES)	for ES cell transfer (for research use only), 20 ° tip angle, 15 µm inner diameter, 1 mm flange, sterile, 25 pcs.	CHF 721.00
5195000087	Piezo Drill Tip ICSI	for piezo-assisted mouse ICSI (for research use only), 25 ° tip angle, 6 µm inner diameter, 6 mm flange	CHF 562.00
5195000095	Piezo Drill Tip ES	for piezo-assisted mouse ES cell transfer (for research use only), 25 ° tip angle, 15 µm inner diameter, 6 mm flange	CHF 562.00

Accessories

Item No.	Description	Price
5196082001	Grip head set 4 size 0, incl.6 o-rings and 2 distance sleeves, for capillaries with outer	CHF 64.68

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 11.05.2026 at 20:53 o'clock.

www.huberlab.ch