

Autoclavable bags, PP LLG-Labware

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Safety & Hygiene](#) / [Disposal](#) / [Disposal bags](#) / [Autoclavable bags, PP LLG-Labware](#)

Product Data

Constructed from very strong PP film with a broad seam at the bottom. For the disposal of contaminated labware. When autoclaving, the neck of the bag must be left open. Film thickness: 50 µm.

Item No.	description	w x l	autoclavable	Price
600.4040.20	highly transparent	200 x 300 mm	at 121 °C	CHF 7.10 CHF 16.10
600.4041.00	---	---	---	CHF 105.00 CHF 234.40
600.4041.01	highly transparent	400 x 780 mm	at 121 °C	CHF 135.00 CHF 300.00
600.4041.02	highly transparent	600 x 800 mm	at 121 °C	CHF 165.00 CHF 370.00
600.4041.03	standard	600 x 800 mm	at 134 °C	CHF 61.00
600.4041.04	highly transparent	700 x 1100 mm	at 121 °C	CHF 175.00 CHF 395.00
600.4041.05	standard	700 x 1100 mm	at 134 °C	CHF 78.90
600.4042.21	standard	400 x 780 mm	at 134 °C	CHF 106.40
600.4042.22	standard	500 x 780 mm	at 134 °C	CHF 128.30
600.4042.23	standard	600 x 780 mm	at 134 °C	CHF 143.90
600.4042.24	standard	700 x 1100 mm	at 134 °C	CHF 254.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 07.05.2026 at 12:21 o'clock.