

# Laboratory Blender Waring® 8010/8011

[Product in eShop](#)

Home / [Labortops\\_Q1\\_2025](#) / Laboratory Blender Waring® 8010/8011

## Product Data

ersatile laboratory mixer ideally suited for general laboratory work. Supplied either with a 1 L stainless steel container or a 1.2 L glass container.

Models 8010 EG/ES with mechanical timer (3 minutes).

- Two velocity settings (18000 or 22000 rpm)
- Epoxy coated motor housing
- Two-piece lid vinyl/styrene

**Scope of supply:** Mixer, container (stainless steel or glass) with handle and lid

- Min. speed: 18000 rpm
- Max. speed: 22000 rpm

Item No.	type no.	description	Price
670.2846.93	8010EG	Glass container, 1.2 l	CHF 840.00 <del>CHF 1'140.00</del>
670.2846.94	8010ES	Stainless steel container, 1 l	CHF 940.00 <del>CHF 1'280.00</del>
670.5714.36	---	---	CHF 935.00 <del>CHF 1'270.00</del>
670.5714.37	---	---	CHF 1'040.00 <del>CHF 1'420.00</del>

## More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 02.05.2026 at 19:24 o'clock.

[www.huberlab.ch](http://www.huberlab.ch)