

Mixer Vortex Genie® - Serie Scientific Industries



Product in eShop

Home / [Labortops_Q1_2025](#) / Mixer Vortex Genie® -Serie Scientific Industries

Product Data

Robust and reliable shakers for speedy mixing and shaking of samples. Multiple accessories: e.g. adapter tops for test tubes, beakers, centrifuge tubes etc. Adapter tops can be changed easily with the snap fit system.

- **Vortex-Genie® 2:** Variable speed for precise mixing from gentle to vigorous. Broad range of attachments for most tubes, plates and other containers. Hands-Free or Touch On control.
- **Vortex-Genie® 2T with an integrated timer:** Timer function: (Touch On) 1-60 seconds, (Hands-Free) 1-60 minutes or continuous. Unit automatically stops at the end of user programmed time.
- **Digital Vortex-Genie® 2 with digital control and display of speed and time:** Digital display provides accuracy, reproducibility and repeatability. Timer function: (Touch On) 1-99 seconds, (Hands-Free) 1-99 minutes or continuous. Specially configured unit for Applied Biosystems protocol available.
- **Vortex Genie® Pulse:** Fully programmable vortex mixer with ON/OFF feature. Pulse 0.1 sec to 99 min.

Scope of supply: Vortexer, platform (diam. 76 mm) with rubber cover, Pop-off cup

Specifications (Vortex-Genie® 2 / Vortex-Genie® 2T / Digital Vortex-Genie® 2 / Vortex Genie® Pulse)

- Shaking movement: orbital
- Speed display/Timer: Scale/No / Scale/Scale / Digital/Digital
- Speed range: 600-2700rpm / 600-2700rpm / 500-2850rpm / 500-2850rpm
- Timer Interval operation: -/ 1 to 60s / 1 to 99s / 0.1 s to 99 min
- Hands-free mode: -/ 1 to 60min / 1 to 99min / 1 to 99 min
- Dimensions (L x W x H): 165 x 122 x 165mm
- Permissible ambient conditions: 0 to 38°C at 95% rh
- Weight: 4kg
- Power supply: 230V, 50/60Hz

17.1378.01	Vortex Genie® 2	CHF 265.00 CHF 485.00
17.1379.10	Vortex Genie® Pulse	CHF 792.45
17.1379.30	Digital Vortex-Genie® 2	CHF 754.65
17.1379.40	Vortex-Genie® 2T	CHF 724.95

Accessories

Item No.	Description	Price
17.1378.04	Accessories for Mixer Vortex Genie® 2 Scientific Industries	from CHF 35.10

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 03.05.2026 at 11:32 o'clock.