



## Pipette tips Low Retention, PP, sterile LLG-Labware

[Product in eShop](#)

Home / [Labortops\\_Q1\\_2025](#) / Pipette tips Low Retention, PP, sterile LLG-Labware

### Product Data

LLG-Pipette tips with hydrophobic inner surface ensure maximum sample recovery and high-precision pipetting. Compatible with all common pipette brands.

- Colour-coded racks facilitate tip identification
- With volume graduation
- Metal-free
- DNase- and RNase-free
- Autoclavable

For research and development only.

- Colour: Transparent
- With graduation: Yes
- Sterile: Yes
- With filter: No
- Low retention: Yes

Item No.	volume	length	tips	Price
600.6944.57	10 µl	31.2 mm	standard	CHF 51.30
600.6944.58	10 µl	45.8 mm	XL	CHF 51.30
600.6944.59	200 µl	53.4 mm	standard	CHF 82.70
600.6944.60	1000 µl	84.3 mm	standard	CHF 42.25

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 22.04.2026 at 18:56 o'clock.