

Serial Extraction Apparatus behrotest® for Soxhlet-/Fat-Extraction Behr

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Environmental technology](#) / [Extraction](#) / Serial Extraction Apparatus behrotest® for Soxhlet-/Fat-Extraction Behr

Product Data

behrotest® serial heaters are the economical and user-friendly solution for the classic Soxhlet or fat extraction:

- complete with reaction vessels, extraction spacers and coolers
- energy for each sample position is individually adjustable
- cooling water distributor ensures even cooling at all sample positions
- practical stand for holding the condensers securely between extractions
- mounts for holding the extraction spacers after removing the sample vessels
- S models: The extractors with stopcock do away with the need for any additional distillation after extraction

Item No.	type no.	description	capacity	Price
624.0430.29	R 104 S	4-position, with stop-cock	100 ml	CHF 3'670.00
624.0430.35	R 304 S	4-position, with stop-cock	30 ml	CHF 3'565.00
624.0430.36	R 306 S	6-position, with stop-cock	30 ml	CHF 4'795.00
624.0430.39	R 604 S	4-position, with stop-cock	60 ml	CHF 3'710.00
624.0430.40	R 606 S	6-position, with stop-cock	60 ml	CHF 5'070.00
624.0430.42	R 254 S	4-position, with stop-cock	250 ml	CHF 4'485.00
624.6023.99	---	---	---	CHF 3'590.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 20:04 oclock.

www.huberlab.ch