

Greenough Stereo Microscopes Lab-Line OZL Kern

[Product in eShop](#)

[Home](#) / [Labortops Q1 2025](#) / [Greenough Stereo Microscopes Lab-Line OZL Kern](#)

Product Data

The flexible all-rounder with zoom function for schools, training companies, inspection authorities and laboratories.

- Tube 45° inclined
- Wide Field HWF 10x/20 mm
- Diopter adjustment: Both-sided
- Illumination: 3 W-LED Incident and 3 W-LED transmitted, separately dimmable
- Interpupillary distance 55 to 75 mm
- Large selection of eyepieces, external illumination units as well as auxiliary objectives are available as accessories
- OZL 466 with strong, continuously dimmable, integrated LED ring illumination in the objective housing
- OZL 468 with integrated handle and stable, mechanical stand
- Protective dust cover, eye piece rubber caps, as well as multi-lingual user instructions included in the scope of supply

Inclination angle tube: 45°

Interpupillary distance: 55 ... 75 mm

Light type: LED

Power supply: 100 ... 240 V

Item No.	type no.	description	enlargement	Price
17.4515.50	OZL 463	Binocular	7x ... 45x	CHF 810.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 22:01 oclock.

www.huberlab.ch