

Sample bags, TWIRL'EM, with writing area, heavy sample closure Labplas

[Product in eShop](#)

Home / [Labortops_Q1_2025](#) / Sample bags, TWIRL'EM, with writing area, heavy sample closure Labplas

Product Data

For collecting, transporting and analysing liquid and solid samples.

- Integrated secure closure system with 1 round wire and 1 flat wire
- Sterile
- Contaminant-free
- Waterproof and airtight
- RNase, DNase and pyrogen free
- With labelling field

Item No.	volume	Dimensions	Thickness	Price
15.1123.08	4000 ml	254 x 382 mm	0.101 mm	CHF 323.75
15.1123.13	450 ml	114 x 229 mm	0.063 mm	CHF 196.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 15:43 o'clock.