

## Sampling spoons, PS, sterile, individually packaged Antonides

[Product in eShop](#)

Home / [Labortops\\_Q1\\_2025](#) / Sampling spoons, PS, sterile, individually packaged Antonides

### Product Data

Designed for precise sampling, this durable polystyrene (PS) spoon is ideal for handling sensitive materials in laboratory or Industrial settings.

- Individually packaged
- Sterilized

Item No.	volume	length	colour	Price
681.7014.02	10 ml	138 mm	white	CHF 37.05
681.7014.03	5 ml	220 mm	red	CHF 67.65

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 21.04.2026 at 10:23 o'clock.