

## Sample containers, PS/PP, with tamper- evident cap, LDPE/PP LLG-Labware



Product in eShop



Home / [Labortops\\_Q1\\_2025](#) / Sample containers, PS/PP, with tamper-evident cap, LDPE/PP LLG-Labware

### Product Data

White. Sturdy containers made of PS or PP. The tamper-evident closure enables a safe transport of samples.

- Food safe
- Not sterile

| Item No.    | volume  | height | diameter         | material       | Price     |
|-------------|---------|--------|------------------|----------------|-----------|
| 600.4023.50 | 40 ml   | 52 mm  | with lid: 36 mm  | PS / lid: LDPE | CHF 14.90 |
| 600.4023.53 | 100 ml  | 58 mm  | with lid: 58 mm  | PP / lid: PP   | CHF 12.60 |
| 600.4023.56 | 250 ml  | 95 mm  | with lid: 64 mm  | PS / lid: LDPE | CHF 17.35 |
| 600.4023.59 | 1000 ml | 116 mm | with lid: 116 mm | PP / lid: PP   | CHF 22.35 |

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 17.04.2026 at 21:58 o'clock.