

Clax folding box, 46 litres Rauschenberger Innovations

[Product in eShop](#)

Home / [Labortops_Q1_2025](#) / Clax folding box, 46 litres Rauschenberger Innovations

Product Data

Sturdy folding box, made of high-quality plastic, for use with the clax Mobil comfort laboratory trolley.

- Weight: 2 kg

Item No.	Volume	Dimensions	colour	Price
633.1480.99	46 l	380 x 540 x 270	grey	CHF 21.00
633.4929.73	46 l	380 x 540 x 270	black	CHF 21.00

Accessories

Item No.	Description	Price
633.1481.00	Laboratory Trolley clax Mobil comfort with Box	from CHF 295.60

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 23:28 oclock.

www.huberlab.ch