

Wireless thermometer XENA HOME TFA



 Product in eShop

Home / [Labortops Q1 2025](#) / Wireless thermometer XENA HOME TFA

Product Data

Wireless thermometer for monitoring the outdoor temperature and indoor climate with temperature and humidity. The outdoor temperature is measured wirelessly via the supplied temperature transmitter. The range is 100 m at a transmission frequency of 433 MHz. Up to 3 temperature transmitters can be linked to the base station.

- Radio-controlled clock with entire date, time zone correction ± 12 hours
- Highest and lowest values are saved with time and date, automatic reset can be set
- For standing or wall mounting

Scope of supply: Base station, temperature transmitter

- Please order batteries separately.
- Int. temperature: 0 ... 50 °C
- Ext. temperature: -40 ... 60 °C
- Dimensions (WxDxH): 82 x 27 x 140 mm
- Weight: 128 g
- Battery/rechargeable battery: 2 x 1.5 V, AA

Item No.	type no.	Price
671.3074.10	XENA HOME	CHF 33.95

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 14:47 o'clock.