

Temperature transmitter with cable sensor TFA



 Product in eShop

Home / [Labortops Q1 2025](#) / Temperature transmitter with cable sensor TFA

Product Data

Transmitter with waterproof cable sensor for measuring and transmitting the temperature wirelessly. The display shows the measured temperature. The range is 100 m at a transmission frequency of 433 MHz.

- Adjustable for channel 1 to 3
- Can be used indoors or outdoors
- For hanging or standing

Batteries are not included.

- Weight: 76 g
- Battery/rechargeable battery: 2 x 1.5 V, AA

Item No.	type no.	description	Dimensions	Price
671.3255.02	Temperature transmitter with cable sensor	Temperature: -40 ... 60 °C	48 x 26 x 129 mm	CHF 13.55

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 15:26 o'clock.