

# Laboratory thermometer LT-105, with glass bottle Dostmann



 Product in eShop

Home / [Labortops\\_Q1\\_2025](#) / Laboratory thermometer LT-105, with glass bottle Dostmann

## Product Data

Ergonomically shaped laboratory thermometer with alarm function for monitoring room temperature, as well as temperature in refrigerators and freezers.

- MAX/MIN function
- Hi-Lo alarm
- Cable probe (external sensor) with approx. 1 m cable
- °C/°F switchable
- With calibration certificate at -18 °C and 4 °C

Laboratory thermometer with cable probe in a 30 ml glass bottle. Incl. Glysofor L sachet for self-filling.

**Scope of supply:** Laboratory thermometer, 1 x AAA battery, holders and holding accessories, 15 ml sachet Glysofor L, calibration certificate, operation manual

- Ext. temperature: -50 ... 70 °C
- Int. temperature: -10 ... 50 °C
- Resolution: 0.1 °C
- Accuracy: ±1 °C
- Sensor (ØxL): 5 x 25 mm
- Cable length: 1 m
- Dimensions (WxDxH): 68 x 40 x 15 mm
- Weight: 43 g

Item No.	type no.	Price
607.6917.94	LT-105	CHF 37.75

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 17.04.2026 at 21:58 oclock.

[www.huberlab.ch](http://www.huberlab.ch)