

Weighing boats, PS, diamond, antistatic LLG-Labware

[Product in eShop](#)

Home / [Labortops_Q1_2025](#) / Weighing boats, PS, diamond, antistatic LLG-Labware

Product Data

Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases. Temperature range from -10 to 70 °C. Can be used with microwave moisture analysers. Features a diamond shape that easily flexes for accurate pouring and a wide, flat bottom to resist tipping.

For research and development only.

Item No.	capacity	colour	Dimensions	Price
600.9007.88	5 ml	white	55 x 35 x 6 mm	CHF 31.65 CHF 89.80
600.9008.01	---	---	---	CHF 34.85 CHF 98.80
600.9501.94	30 ml	white	85 x 65 x 14 mm	CHF 12.40 CHF 35.80
600.9501.95	100 ml	white	125 x 100 x 19 mm	CHF 18.10 CHF 44.80
600.9501.96	---	---	---	CHF 13.70 CHF 39.30
600.9501.97	---	---	---	CHF 19.35 CHF 49.30

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 03.05.2026 at 11:31 o'clock.