

Docking station DST 01 for Ultrasonic baths SONOREX smart ST Bandelin



Product in eShop

Home / [Labortops Q2 2025](#) / Docking station DST 01 for Ultrasonic baths SONOREX smart ST Bandelin

Product Data

Easy transfer of the process protocol to other storage media.

- USB-A, USB-B and Ethernet connection
- IP code: IP 20

Item No.	type no.	Price
17.9967.33	DST 01 for SONOREX smart ST	CHF 314.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 03.05.2026 at 11:31 o'clock.