

Photometer MD50 COD Tintometer Lovibond

NEW



 Product in eShop

Home / [Labortops_Q3_2025](#) / Photometer MD50 COD Tintometer Lovibond

Product Data

The MD50 COD is specifically designed for COD determination and offers superior performance, enhanced usability and expanded testing capabilities, including the new Very Low Range (VLR) COD method for ultra-sensitive measurements.

- Highest accuracy for low-pollution samples
- Measurement of COD in high, medium and very low ranges for versatile applications
- Easy operation thanks to intuitive symbols and animations as well as user-friendly navigation
- Results within seconds enable high efficiency
- Time saving with the One Time Zero function (stays valid for up to 8 hours)
- Handheld design for comfortable single-handed operation
- Robust housing with side protection for extra durability
- Wireless NFC and USB-C data transfer
- Automatic GLP compliance
- Versatile connectivity via the AquaLX® app or Lovibond® Data Expert software
- Pre-programmed Lovibond® CSB cuvette tests for operation right out of the box available separately

Methods: M134, M130, M133, M131, M132

Scope of supply: Photometer in carrying case, adapter for cells (16 mm), 2 x dosing syringes (1 ml, 2 ml), USB-C cable, screwdriver, manual, 2 x AA batteries

- Photometric measuring range: 0 ... 3 Abs
- Illuminants: LED
- Dimensions (WxDxH): 76 x 155 x 45 mm
- Weight: 247 g
- Battery/rechargeable battery: 2 x 3 V, AA, LR6, AIMn
- IP code: IP 67

Item No.	type no.	Price
675.5236.19	MD50 COD	CHF 380.00
		CHF 495.00

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 26.04.2026 at 15:16 o'clock.

www.huberlab.ch