

## Bowl Trolleys GORR



Product in eShop

Home / [Labortops Q3 2025](#) / Bowl Trolleys GORR

### Product Data

#### Trolley

- Blue-grey solid rubber (non-marking) tyres
- Welded frame for 85 or 100 litres plastic bowls
- Height with red bowl approx. 760 mm, with white bowl approx. 840mm
- 4 swivel castors with 75 mm diam.
- ball bearing mounted
- powder coated
- Carrying weight approx. 100 kg

**Bowl, red:** PE, 85 litres**Bowl, white:** PE, 100 litres

| Item No.    | description                 | colour       | Dimensions      | Price      |
|-------------|-----------------------------|--------------|-----------------|------------|
| 692.0842.01 | Tank trolley (without tank) | black-silver | 500 x 750 x 720 | CHF 271.70 |
| 692.0842.02 | Tank                        | red          | 540 x 765 x 350 | CHF 121.30 |
| 692.0842.03 | Tank                        | white        | 540 x 900 x 400 | CHF 153.75 |

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 21.04.2026 at 22:01 o'clock.