

## Support rods with thread galvaniser steel Bochem

[Product in eShop](#)

Home / [Labortops\\_Q3\\_2025](#) / Support rods with thread galvaniser steel Bochem

### Product Data

Steel, galvanised with thread

- M10 thread
- Material: Steel
- Diameter: 12 mm

Item No.	length	Price
13.5150.75	500 mm	CHF 19.70
13.5151.76	600 mm	CHF 21.20
13.5152.77	750 mm	CHF 23.65
13.5153.78	1000 mm	CHF 26.90
13.5154.79	1250 mm	CHF 30.30
13.5155.80	1500 mm	CHF 34.40

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 22.04.2026 at 11:11 o'clock.