

Test tube shaker LLG-uniTEXER 1 pro LLG-Labware

[Product in eShop](#)

Home / [Labortops_Q3_2025](#) / Test tube shaker LLG-uniTEXER 1 pro LLG-Labware

Product Data

Compact test tube shaker with outstanding price-performance ratio and new, improved design. Ideal for mixing small sample volumes in vessels up to 30 mm diameter.

- Touch function for short-time operation
- Adjustable speed control for excellent mixing results
- Suction cups ensure a secure hold during vortexing
- Extremely small footprint

Scope of supply: Test tube shaker, universal power supply incl. adapter for EU, UK, US, AUS

- Type of movement: orbital
- Shaking stroke: 4.5 mm
- Speed: 0 ... 4,200 1/min
- Dimensions (ØxH): 120 x 103 mm
- Weight: 750.0 g
- Permitted ambient temperature: 4 ... 45 °C
- Rated capacity: 9 W
- Plug type: Universal
- Power supply: 100 ... 240 V, 50/60 Hz
- IP code: IP 43
- Warranty: 3 years

Item No.	type no.	Price
600.2635.09	LLG-uniTEXER 1 pro	CHF 166.65

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 18:23 oclock.

www.huberlab.ch