

# Conductivity Total Dissolved Solids TDS Standards Reagecon

[Product in eShop](#)

Home / [Reagecon\\_ConductivityStandards](#) / Conductivity Total Dissolved Solids TDS Standards Reagecon

## Product Data

The Total Dissolved Solids (TDS) Standards are aqueous based, thereby eliminating any errors attributable to matrix mismatch. Each bottle is supplied with a comprehensive Certificate of Analysis. The low temperature coefficient of variation, reduces measurement errors and enables non-temperature controlled use and is suitable for field work. Customised or bespoke values can be manufactured on demand.

Item No.	Content	Packaging	Price
500.CS1382	500 milliliter	HDPE Bottle GL45	CHF 51.60
500.TDS10	1 liter	HDPE Bottle GL45	CHF 125.10
500.TDS100	1 liter	HDPE Bottle GL45	CHF 83.30
500.TDS1000	1 liter	HDPE Bottle GL45	CHF 59.40
500.TDS1500	1 liter	HDPE Bottle GL45	CHF 59.40
500.TDS2000	1 liter	HDPE Bottle GL45	CHF 59.40
500.TDS250	1 liter	HDPE Bottle GL45	CHF 141.65
500.TDS50	1 liter	HDPE Bottle GL45	CHF 141.65
500.TDS500	1 liter	HDPE Bottle GL45	CHF 141.65

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 20.04.2026 at 14:47 o'clock.