

Buffer Standard pH 4.00/7.00 at 20/25°C Colour Coded ISO 17034 Certified Reference Material Reagecon



A CALIBRE SCIENTIFIC COMPANY



Product in eShop

Home / Brands / Reagecon Diagnostics / Buffer Standard pH 4.00/7.00 at 20/25°C Colour Coded ISO 17034 Certified Reference Material Reagecon

Product Data

Reagecon's Buffer Standard Certified Reference Material is colour coded for ease of use, produced under our ISO 17034 (A2LA Ref: 6739.01) Accreditation and is NIST traceable.

Item No.	Content	Packaging	Price
500.CRM1040525C	1 piece		CHF 23.00 CHF 32.60
500.CRM10405C	1 piece		CHF 22.50 CHF 31.60
500.CRM1070525C	1 piece		CHF 23.00 CHF 32.60
500.CRM10705C	1 piece		CHF 22.50 CHF 31.60

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 12.05.2026 at 07:20 o'clock.