

Blood Collection Tray

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Medicine & Healthcare](#) / [Hematology](#) / [Blood Collection](#) / Blood Collection Tray

Product Data

Light weight, heavy on organization!

Organize your phlebotomy, I.V. start or catheter insertion supplies with this spacious, robust and easy to clean carry tray. Customize the compartments to suit your personal needs with the dividers provided.

- Lightweight easy-to-clean ABS construction
- Comfort grip handle tested to hold up to 10 Kg!
- Comes complete with a 17 mm well 60 place tube rack
- Tray Measures 419 x 311 x 85 mm, handle height 241 mm
- Tube rack 17 mm

Item No.	Price
15.1202.74	CHF 172.55

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 05:43 o'clock.