

Microscope slide box

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Medicine & Healthcare](#) / [Microscopy](#) / [Slides](#) / Microscope slide box

Product Data

Made of ABS, for 100 microscope slides of size 76 x 52 mm. With hinged lid, nickel-plated catch and cork inlay. Index card integrated in the lid. The numbering of the index card matches the numbering embossed on the base of the box.

- Ideal for storage and transport
- Stackable

Item No.	l x w x h	Price
15.1597.39	208 x 175 x 60 mm	CHF 38.35

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 18.04.2026 at 22:23 o'clock.