

Embedding cassettes, Histosette® II Standard

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Microscopie](#) / Embedding cassettes, Histosette® II Standard

Product Data

Made of acetal polymer

- Resistant to the chemical action of histological solvents
- Flow-through slots maximise fluid exchange
- Integral lid
- Anterior writing area at 45° angle to enable use with leading cassette labelling instruments

Item No.	colour	Price
10.9492.02	white	CHF 250.35
10.9492.03	pink	CHF 250.35
10.9492.04	green	CHF 250.35
10.9492.05	yellow	CHF 250.35
10.9492.06	blue	CHF 250.35
10.9492.07	peach	CHF 250.35
10.9492.08	beige	CHF 250.35
10.9492.09	gray	CHF 250.35
10.9492.10	mauve	CHF 250.35
10.9492.11	orange	CHF 250.35
10.9492.12	aqua	CHF 250.35

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 07.05.2026 at 12:23 oclock.

www.huberlab.ch