

## Deep well plate, 384-well, stackable 300 µl, PP, non-sterile, pack of 48 alphanumeric

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Life Science](#) / [Highthroughput Screening \(HTS\)](#) / [Microplates](#) /

### Product Data

Aus PP, Volumenbereich 0,3 ml. Für eine Vielzahl von Anwendungen entwickelt, z.B. High-Throughput Screening (HTS), allgemeine Untersuchungen, die eine Stammplatte benötigen, Zell- und Gewebekultur, Reihenverdünnungen, Reagenztransfer, Probenaufbewahrung bis zu -80°C etc. Deep-well Platten sind im ANSI/SLAS-Format hergestellt und können platz sparend gestapelt werden.

Item No.	Price
7.7013.55	CHF 585.95

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 07.05.2026 at 03:32 o'clock.