

Modular glove box holders



Product in eShop

[Home](#) / [Assortment](#) / [Medicine & Healthcare](#) / [Modular glove box holders](#)

Product Data

Modular design allows the mounting of a single bracket or the connected assembly on top of each other.

Colour blue: with Microban® coating that inhibits the growth of mould and bacteria. Ideal for use in clean rooms, hospitals and microbiological work areas.

- Universal fit for the most common glove boxes
- The highly impact-resistant, durable resin is resistant to harsh cleaning solutions
- Modular design allows single boxes or horizontal stacking for multiple box holders
- Screws for mounting are included in the delivery
- Dimensions: 254 x 156 x 97 mm

Item No.	colour	Price
3.7899.01	white	CHF 108.45
3.7899.02	blue	CHF 177.25

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

www.huberlab.ch