

Mobile Overshoe Dispenser HYGOSTAR Franz Mensch



Product in eShop

Home / [Labortops_Q3_2025](#) / Mobile Overshoe Dispenser HYGOSTAR Franz Mensch

Product Data

Overshoe dispenser HYGOSTAR

For employees and visitors to hotels, companies, schools, laboratories, offices, catering facilities, hospitals, dairies, building trades, exhibition rooms and buildings, production facilities, clean rooms and protection zones.

Overshoe dispenser Easy Box

- Particularly cost-effective overshoe dispenser
- for low visitor numbers
- Space for 40 overshoes
- Very space-saving
- Flexible use, as no electricity required

Overshoe dispenser HYGOMAT® Classic

No bending down, with handle to hold on to - for even more comfortable donning, four times faster than normal overshoe donning, with viewing window so that the fill level can be seen at all times, refilled with overshoes in just a few seconds.

- 4x faster than putting on overshoes manually
- Space-saving and without electricity
- Space for up to 110 overshoes
- easy handling, very quiet

Overshoe dispenser HYGOMAT® Comfort

No bending down, with handle to hold on to - for even more comfortable donning, four times faster than normal overshoe donning, with viewing window so that the fill level can be seen at all times, refilled with overshoes in just a few seconds.

- Extremely efficient system for disposable overshoes
- 4x faster than manual donning
- very user-friendly, no bending down
- works without electricity
- With handle for easy donning
- Space for up to 220 overshoes

Overshoe dispenser HYGOMAT® Cleanroom

Particularly hygienic as the back and side panels are made of stainless steel, making it very easy to clean and disinfect, with 2 castors at the back - can be quickly and easily rolled to another location, very flexible, especially for cleanrooms and pharmaceuticals.

- 4x faster than putting on overshoes manually
- Suitable for ISO class 5 cleanrooms
- Particularly hygienic as it is made of stainless steel (back and side panels)
- can therefore be cleaned and disinfected very easily
- Flexible application options, as no electricity is required
- extremely clean, as overshoes are not contaminated by touching them with your hands
- 100% hygienic, no cross-contamination

Item No.	description	Dimensions	weight	Price
3.4924.06	Overshoe dispenser Easy Box	400 x 210 x 120 mm	1.5 kg	CHF 84.00
3.4924.08	Overshoe dispenser HYGOMAT® Classic	695 x 315 x 230 mm	8.5 kg	CHF 366.00
3.4924.09	Overshoe dispenser HYGOMAT® Cleanroom	740 x 300 x 730 mm	17.5 kg	CHF 725.00
3.4924.11	Overshoe dispenser HYGOMAT® Comfort	740 x 300 x 730 mm	17,5 kg	CHF 587.00

Accessories

Item No.	Description	Price
3.4924.13	Overshoes for Overshoe Dispenser HYGOMAT® Franz Mensch	from CHF 195.58

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 02.05.2026 at 06:25 oclock.

www.huberlab.ch