

# Syringe SteriPlast

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Liquid Handling](#) / [Sampling](#) / [Syringe SteriPlast](#)

## Product Data

The syringe is used to take sterile samples and store them. The sample amount can be read off from the printed contents graduation. The syringe is supplied in a sterile package and is ready for immediate use. The syringe can be extended to any length with additional tubing. A Luer adapter is supplied for the 100 ml syringe.

- PP, transparent
- Contents graduation
- Sterilised, packed separately

Item No.	vol.	division	Price
15.5325.06	50 ml	1 ml	CHF 29.00
15.5325.10	100 ml	2 ml	CHF 86.20

## More product images





\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 09.05.2026 at 15:49 o'clock.

[www.huberlab.ch](http://www.huberlab.ch)